

Accessories for M2 shower/commode chair Increased and reduced armrests



Armrest with increased width

If you want an increased width between the armrests it can be achieved by acquiring a set of armrest.

One armrest gives an increase of 3 cm.

Armrest with reduced width

If you need the armrests to be closer to the body you can purchase a set of armrests that makes the width 5 cm smaller.

One armrest reduce the space with -2.5 cm



HMN since 1934

Technical data

Increased armrests

Armrest with increased width

Fits on all M2 chairs

Article no. 1 set increased armrests +6 cm	310390
Article no. For right armrest + 3 cm	310391
Article no. For left armrest + 3 cm	310392

Reduced armrests

Armrest with reduced width

Fits on all M2 chairs

Article no. 1 set reduced armrests -5cm	310393
Article no. For right armrest – 2.5 cm	310394
Article no. For left armrest – 2.5 cm	310395



HMN since 1934